



IMPORTANT SAFETY INSTRUCTIONS

Read all instructions carefully before use !

- § Risk of electric shock, may interfere with electrical devices such as pacemakers
- § This high voltage instrument is not intended for clinical measurements using human subjects.
- § Do not touch the **Output** connector either directly or indirectly e.g. through a wire connected to the connector while the device is **on**. Before any manipulation with the wires connected to the **Output** connector make sure that the device is turned **off** – the device switch must be in the **OFF** position and the blue light under the switch must be off.
- § Do not open the device, replace the internal battery, or connect different accessories than the ones specified in the device documentation or use any other charger than the

one provided by Bio-Signal Group Corp. In case of problems operating the device, contact Bio-Signal Group Corp for repairs. Do not attempt to repair the device yourself.

- § As with all electrical devices, use caution when operating the device. Make sure the device is turned off before you plug or unplug the device. Keep the device away from liquids
- § Bio-Signal Group Corp. does not assume responsibility for injury or damage due to the misuse of this instrument.

